
MPI

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Mandibular Position Indicator

(MPI) is used to register the three dimensional position of the mandible in relation to the cranium.

Occlusal interferences and / or disturbances of the neuromuscular System often cause changes in condylar position.

These changes normally occur between retruded contact position (RCP) and intercuspal position (IP) or the habitual position.


The magnitude and direction of these changes can be measured with the MPI. The Mandibular Position Indicator (MPI) consists of a modified upper member of the SAM® articulator with sliding cubes instead of condylar housings. The measurements are made between the condylar elements of the lower member and the sliding cubes of the upper member.

The patented MPI records the measurements on the transverse axis within the intercondylar distance of the condylar elements of the articulator. Therefore, this system has eliminated reciprocal readings and misinterpretations that are characteristic of measurements made outside of condylar elements.

Since the three dimensional measurements are all made to the same reference plane, you are able to make direct interpretations of the data. This data can be directly compared before, during, or after therapy or at any other time. The data can be transferred directly to the cephalometric xray so that articulator information can be superimposed on an x-ray.

Fig. 1

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
**Fig. 1**

shows incisal pin deviation due to reciprocal horizontal rotation in a typical lateral movement (top view).

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
**Fig. 2**

shows the same incisal pin deviation due to pure translation (side shift) of the lower member of the articulator (top view).

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**Fig. 3**

shows the same incisal pin deviation due to a pendulum effect in a reciprocal vertical rotation of the condylar elements (frontal view).

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Procedure:

1. *Mount the maxillary Cast with the SAM® Anatomical Transfer Bow or Hinge Axis Transfer Bow (SAM AXIOGRAPH®) and then mount the mandibular cast with a Centric Inter-occlusal Record. Verify mounting.*
2. *Place the incisal table on the lower member of the articulator and the incisal pin assembly on the upper member.*
3. *Lock the articulator in RCP with the Centric Locking Devices and raise the incisal pin.*
4. *Attach the incisal graph paper to the incisal table and place RED marking tape on top of the graph paper with the marking side towards the paper.*
5. *Close the articulator carefully in terminal hinge relationship (RCP) and move the incisal pin to lightly come in contact with the tape and graph paper. Tap the pin one time to get a clean mark on the graph paper. Lock the incisal pin.*
6. *Read the height value of the incisal pin and enter the value on the MPI record form under RCP.*
7. *Remove the upper member of the articulator and transfer the maxillary cast to the Mandibular Position Indicator (MPI), (Fig.4).*
8. *Prepare the MPI by inserting the incisal pin halfway (Fig.5).*

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Release the gauge arm and turn the arm to an upright position. Hold the cubes away from the frame to avoid perforating the graph paper prematurely. Place the condylar graphs parallel to the edges of the cubes and use the articulator symbols on the graph as your guide for right and left.

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Do not perforate the graph papers at this time. cap is on top of pin left. Do not perforate the graph papers at this time.

NOTE

With the SAM 2 & 3 Incisal Table, the red rings on the incisal pin are up above the upper member; above the millimeter scale. Black plastic cap is on top of pin above red rings.

With the SAM 1 Incisal Table, the red rings on the incisal pin are down toward the incisal

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9. *Carefully position the maxillary cast into the IP position and hold securely against the mandibular cast (Fig.6).*

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
10. *Place BLACK marking tape on top of the graph paper on the incisal table. Care-fully lower the incisal pin, tap one time against the incisal pin and lock the incisal pin. Read the height value of the incisal pin and enter the value on the MPI record form under IP.*

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11. *Hold the MPI assembly firmly in position and place BLACK marking tape between the sliding cube and condylar element with the marking side towards the cube. Move the cube and tap against the condylar element to obtain a recording of the IP position in the horizontal (X) and vertical (Z) position. Repeat the procedure for the opposite side.*

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
The transverse displacement of the mandible is measured and is identified as the Y axis and is the lateral movement of the black cubes.

There is 5 mm of space between the inner wall of the black cubes and the MPI frame.


Since the distance is equal on both sides, it is only necessary to measure one side and determine whether it moves medially or laterally.

The left side is measured and identified with positive for lateral movement and negative for medial movement.

Swing the gauge arm down into the left sliding cube and insert the arm to the depth of the slot. The dial gauge registers the clearance between the MPI frame and the left sliding cube. Increased space indicates positive readings, black numbers on the dial gauge, and the mandible is shifted to the left. Decreased space indicates

$$\begin{aligned} \text{Delta I} &= 3 \\ \text{Delta y} &= 0,4 \\ \text{Delta z} &= \underline{0 \mid +3} \\ \text{Delta x} &= \underline{+2 \mid 0} \end{aligned}$$

$$\begin{aligned} \text{IP} &= 0 \\ \text{RCP} &= -2 \\ \text{Delta h} &= \underline{-2} \end{aligned}$$

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negative readings, red numbers on the dial gauge, and the mandible is shifted to the right (Fig. 7).

The change between RCP and habitual occlusion position can also be measured by inserting a habitual occlusion record between the casts and using a different colored tape. Repeat the measurements. Disengage the dial gauge arm, swing it upward away from the MPI assembly, and slide the black cube inward. A needle will penetrate the recording graph to represent the condylar mounting position of the centric interocclusal record. Repeat the procedure for the opposite side.

NOTE

Perforation prior to this point will distort the recording surface and give poor recordings. A ring or semilunar mark around the perforation indicates coincidence with the centric mounting position.

Adjustment:

Swing the gauge arm down into the left sliding cube and insert the arm to the length of the slot. Insert the 50x5mm cylinder between the frame and the left cube below the axis. Hold the sliding cube firmly against the cylinder. The pointer on the gauge should be at the zero mark.

If necessary, reset the hexagon set-screw in the cube using the enclosed allen wrench.

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*Mandibulator Position Indicator Labels
(package) 15 sheets with 6 sets per sheet.
For use with MPI.*

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Mandibulator- Position- Indicator, with optical loupe used for comparing centric records, changes in mandibular cast position, monitoring splint adjustments, etc. Includes universal incisal pin, labels, calibration tool and optical loupe.

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Like MPI 316M, for systems with screw type mounting plates

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Condymeter III, analog MPS

3-dimensional dial gauge version of MPI, accuracy 0.01mm, primarily for research and teaching, used with screw type mounting plate or AXIOSPLIT® System.

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Condymeter III, analog like MPI 415M, for systems with screw type mounting plates

US PAT 4.459.108
e brevetti in altri paesi / *and PAT
in other countries*

Letteratura / Literature

McCollum, B.: Fundamentals involved in prescribing restorative remedies. Dent. Items 61, 522, 641, 724, 852, 942 (1939).

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Hoffmann, P., et al.: Comparison of condylar position in centric relation and centric occlusion. J. prosth. Dent. 30, 582-588, (1973).

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Long, J.: Location by the Terminal Hinge Axis by intraoral means. J. prosth. Dent. 23,11-24 (1970).

Lundeen, H.: Centric relation records. J. prosth. Dent. 31,244-253(1974).

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Sears, V.: Mandibular condyle migrations as influenced by tooth occlusion. J. Amer. dent. Ass. 45, 179-192 (1952).

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Per la formazione pratica in prima persona
For lectures or practical courses in our advanced education institut

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I nostri prodotti vengono fabbricati in Germania.
Our products are all Made in Germany.

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info@sam-dental.de • www.sam-dental.de

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The logo features a large, stylized blue letter 'D' above the words 'DENTAURUM' and 'ITALIA' in a bold, sans-serif font. A yellow horizontal bar is positioned above the logo.

Dentaurum Italia S.p.a. | Via degli Speciali, 142/144 | Centergross 40050 Funo (BO)
Telefono 051 862580 | Fax 051 863291 | www.dentaurum.it | info@dentaurum.it

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are registered trademarks of SAM® Präzisionstechnik GmbH
sono marchi registrati di SAM® Präzisionstechnik GmbH.
SAM® Präzisionstechnik GmbH