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Bibliografia

1. Chen ST, Buser D. Clinical and esthetic outcomes of implants placed in postextraction sites. *Int J Oral Maxillofac Implants*. 2009;24(Suppl):186-217.
2. Lang NP, Pun L, Lau KY, Li KY, Wong MC. A systematic review on survival and success rates of implants placed immediately into fresh extraction sockets. *Clin Oral Implants Res*. 2012;23(Suppl 5):39-66.
3. Buser D, Chappuis V, Belser UC, Chen S. Implant placement post extraction in esthetic single tooth sites: when immediate, when early, when late? *Periodontol 2000*. 2017;73(1):84-102.
4. Chappuis V, Araújo MG, Buser D. Clinical relevance of dimensional bone and soft tissue alterations post-extraction in esthetic sites. *Periodontol 2000*. 2017;73(1):73-83.
5. Kan JY, Rungcharassaeng K, Lozada JL, Zimmerman G. Facial gingival tissue stability following immediate placement and provisionalization of maxillary anterior single implants. *J Prosthet Dent*. 2011;106(6):371-377.
6. Esposito M, Grusovin MG, Polyzos IP, Felice P, Worthington HV. Timing of implant placement after tooth extraction: immediate, immediate-delayed or delayed implants? *Cochrane Database Syst Rev*. 2010;(3):CD005968.
7. De Rouck T, Collys K, Cosyn J. Immediate single-tooth implants in the anterior maxilla: a 1-year prospective study on hard and soft tissue response. *J Clin Periodontol*. 2008;35(7):649-657.
8. Cosyn J, Eghbali A, De Bruyn H, Dierens M, De Rouck T. Single immediate implants in the anterior maxilla: 3-year results of a case series on hard and soft tissue response and aesthetics. *J Clin Periodontol*. 2011;38(8):746-753.
9. Botticelli D, Berglundh T, Lindhe J. Hard-tissue alterations following immediate implant placement in extraction sites. *J Clin Periodontol*. 2004;31(10):820-828.
10. Sanz M, Cecchinato D, Ferrus J, Pjetursson EB, Lang NP, Lindhe J. A prospective, randomized-controlled clinical trial on implants placed in fresh extraction sockets: soft tissue outcomes. *Clin Oral Implants Res*. 2014;25(2):185-193.
11. Morton D, Gallucci G, Lin WS, et al. Group 2 ITI Consensus Report: Prosthodontics and implant dentistry. *Clin Oral Implants Res*. 2018;29(Suppl 16):215-223.
12. Urban IA, Monje A, Lozada JL, Wang HL. Long-term evaluation of peri-implant bone level after immediate implant placement with guided bone regeneration. *Int J Oral Maxillofac Implants*. 2017;32(1):137-144.