




Gli assortimenti

Brackets linguali

Caso	REF	Contenuto	Brackets*	REF
1. Confezione Social six	891-406-00	6 x	Brackets anteriori	790-450-00
2. Arcata singola	891-407-00	6 x	Brackets anteriori	790-450-00
		4 x	Brackets posteriori	790-452-00
3. Arcata singola + tubi	891-408-00	6 x	Brackets anteriori	790-450-00
		4 x	Brackets posteriori	790-452-00
		2 x	Tubi molari	724-102-51
4. Caso completo	891-409-00	12 x	Brackets anteriori	790-450-00
		8 x	Brackets posteriori	790-452-00
		4 x	Tubi molari	724-102-51

* Slot: $\varnothing$ 0,51 x 0,51 mm / 20 x 20



La sequenza degli archi

Archi linguali

Caso	No.	Diametro	Materiale	REF
1. Social six	1.	ø 0,38 mm / 15	dentaflex® archi linguali, tondi, 6x intrecciati	766-771-00
	2.	ø 0,25 mm / 10	rematitan® "LITE" archi linguali, tondi	766-760-00
	3.	ø 0,30 mm / 12	rematitan® "LITE" archi linguali, tondi	766-761-00
	4.	ø 0,40 mm / 16	remanium® archi linguali, tondi	766-767-00
2. Arcata singola	1.	ø 0,38 mm / 15	dentaflex® archi linguali, tondi, 6x intrecciati	766-771-00
	2.	ø 0,25 mm / 10	rematitan® "LITE" archi linguali, tondi	766-760-00
	3.	ø 0,30 mm / 12	rematitan® "LITE" archi linguali, tondi	766-761-00
	4.	ø 0,35 mm / 14	rematitan® "LITE" archi linguali, tondi	766-762-00
	5.	ø 0,40 mm / 16	rematitan® "LITE" archi linguali, tondi	766-763-00
	6.	ø 0,40 mm / 16	remanium® archi linguali, tondi	766-767-00
3. Arcata singola + tubi	1.	ø 0,38 mm / 15	dentaflex® archi linguali, tondi, 6x intrecciati	766-771-00
	2.	ø 0,25 mm / 10	rematitan® "LITE" archi linguali, tondi	766-760-00
	3.	ø 0,30 mm / 12	rematitan® "LITE" archi linguali, tondi	766-761-00
	4.	ø 0,35 mm / 14	rematitan® "LITE" archi linguali, tondi	766-762-00
	5.	ø 0,40 mm / 16	rematitan® "LITE" archi linguali, tondi	766-763-00
	6.	∅ 0,41 x 0,41 mm / 16 x 16	rematitan® "LITE" archi linguali, quadrati	766-765-00
	7.	∅ 0,46 x 0,46 mm / 18 x 18	remanium® archi linguali, quadrati	766-769-00
4. Caso completo	1.	ø 0,38 mm / 15	dentaflex® archi linguali, tondi, 6x intrecciati	766-771-00
	2.	ø 0,25 mm / 10	rematitan® "LITE" archi linguali, tondi	766-760-00
	3.	ø 0,30 mm / 12	rematitan® "LITE" archi linguali, tondi	766-761-00
	4.	ø 0,35 mm / 14	rematitan® "LITE" archi linguali, tondi	766-762-00
	5.	ø 0,40 mm / 16	rematitan® "LITE" archi linguali, tondi	766-763-00
	6.	∅ 0,41 x 0,41 mm / 16 x 16	rematitan® "LITE" archi linguali, quadrati	766-765-00
	7.	∅ 0,48 x 0,48 mm / 19 x 19	rematitan® "LITE" archi linguali, quadrati	766-766-00
	8.	∅ 0,48 x 0,48 mm / 19 x 19	remanium® archi linguali, quadrati	766-770-00